



## Mouse anti-Actin, smooth muscle

Cat. No.: BMS001 (16 ml Ready-to-use)

### Instructions for use

#### Intended use

This antibody is designed for the specific localisation of Actin (smooth muscle actin) in formalin-fixed, paraffin-embedded tissue sections.

Anti-Actin antibody is intended for in vitro diagnostic use.

#### Specifications

<b>Specificity:</b>	Actin (smooth muscle actin)
<b>Immunogen:</b>	Recombinant peptide according to the 10 N-terminal amino acid residues of smooth muscle actin acetylated at the N-terminus
<b>Clone:</b>	1A4
<b>Isotype:</b>	Mouse IgG2a kappa
<b>Species reactivity:</b>	Human +, mouse +, rat +, others not tested

#### Summary and Description

Anti-Actin antibody stains positive in vascular smooth muscles as well as smooth muscle cells of the gastrointestinal tract and the myometrium. Additionally, myoepithelial cells of breast and sweat glands are stained.

The antibody detects tumours derived from smooth muscles and myoepithelial cells.

#### Reagent provided

Mouse monoclonal antibody in buffer with carrier protein and preservative for stabilisation in the following formats:

**Ready-to-use:** 16 ml (Cat. No. BMS001)

#### Dilution of primary antibody

None

#### Storage and handling

The antibody should be stored at 2-8°C without further dilution.

If necessary, dilutions of the antibody should be done with a suitable antibody dilution buffer (e.g. ZUC025 from Zytomed Systems). The diluted antibody should be stored at 2-8°C after use. Stability of this working solution depends on various parameters and has to be confirmed by appropriate controls.

The antibody provided is suitable for use until the expiry date indicated on the label, if stored at 2-8°C. Do not use product after the expiry date. Positive and negative controls should be run simultaneously with all specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Zytomed Systems' technical support or your local distributor.

#### Precautions

Use through qualified personnel only. Wear protective clothing to avoid contact of reagents and specimens with eye, skin and mucous membranes. If reagents or specimens come in contact with sensitive area, wash with large amounts of water.

Microbial contamination of the reagent must be avoided, since otherwise non-specific staining may occur.

ProClin300 and sodium azide (NaN<sub>3</sub>) are used for stabilisation. Reaction of sodium azide with lead or copper in drainage pipes can result in the formation of highly explosive metallic azides. Discard the antibody solution in a large volume of running water to avoid formation of deposits. A material safety data sheet (MSDS) for the pure substances is available upon request.

## Staining procedure

Refer to the following table for conditions specifically recommended for this antibody. Also refer to detection system data sheets for guidance on specific staining protocols or other requirements.

<u>Parameters</u>	<u>Zytomed Systems recommendations</u>
*Pre-treatment	None
*Control tissue	Smooth muscle or leiomyosarcoma
*Working dilution	None
*Incubation time	30 - 60 minutes

## Quality control

The recommended positive control tissues for this antibody are smooth muscles or a leiomyosarcoma. We recommend carrying out a positive and a negative control with every staining run. Please refer to the instructions of the detection system for guidance on general quality control procedures.

## Troubleshooting

If you observe unusual staining or other deviations from the expected results please read these instructions carefully, refer to the instructions of the detection system for relevant information or contact your local distributor.

## Expected results

This antibody stains positive in the cytoplasm of smooth muscle cells in formalin-fixed, paraffin-embedded tissue sections. Interpretation of the staining results is solely the responsibility of the user. Any experimental result should be confirmed by a medically established diagnostic procedure.

## Limitations of the Procedure

Immunohistochemistry is a complex technique involving both histological and immunological detection methods. Tissue processing and handling prior to immunostaining, for example variations in fixation and embedding or the inherent nature of the tissue can cause inconsistent results (Nadji and Morales, 1983). Endogenous peroxidase, alkaline phosphatase or biotin may cause non-specific staining depending on the detection system used. Tissues containing Hepatitis B Surface Antigen (HBsAg) may give false positive results with HRP (horse radish peroxidase) detection systems (Omata *et al*, 1980). Inadequate counterstaining and mounting can influence the interpretation of the results.

Zytomed Systems warrants that the product will meet all requirements described from its shipping date until the expiry date is reached, if the product is stored and utilised as recommended. No additional guarantees can be given. Under no circumstances shall Zytomed System be liable for any damages arising out of the use of the reagent provided.

## Performance characteristics

Zytomed Systems has conducted studies to evaluate the performance of the antibody for use with a standard detection system. The product has been found to be sensitive and specific to the antigen of interest with minimal or no cross-reactivity.

## Bibliography

Alvi AJ *et al*. Breast Cancer Res 5:R1-R8, 2003  
Ribeiro-Silva A *et al*. Arch Pathol Lab Med 127:336-340, 2002  
Reis.Filho JS, Schmitt FC. Adv Anatom Pathol 9:280-289, 2002






Nadji M and Morales AR Ann N.Y. Acad Sci 420:134-9, 1983  
Omata M *et al*. Am J Clin Pathol 73(5): 626-32, 1980

October 22, 2015

Rev: 1015

Doc: DBE\_BMS001

Explanations of the symbols on the product label:

<b>REF</b>	Bestellnummer Catalog Number Reference du catalogue		Verwendbar bis Use By Utiliser jusque		Gebrauchsanweisung beachten Consult Instructions for use Consulter les instructions d'utilisation
<b>LOT</b>	Chargenbezeichnung Batch Code Code du lot		Lagerungstemperatur Temperature Limitation Limites de température	<b>RUO</b>	Nur für Forschungszwecke For Research Use Only Pour la recherche uniquement
<b>IVD</b>	In vitro Diagnostikum In Vitro Diagnostic Medical Device Dispositif médical de diagnostic in vitro		Achtung Warning Attention		Hersteller / Manufacturer / Fabricant Zytomed Systems GmbH • Anhaltinerstraße 16 14163 Berlin, Germany • Tel: (+49) 30-804 984 990 www.zytomed-systems.com