



ZytoChem Plus Double Stain Polymer Kit

REF / Cat. No.: **POLDS-006** **60 tests, 6 ml**
POLDS-100 **1000 tests, 100 ml**

Instructions for use

Intended use

ZytoChem Plus Double Stain Polymer Kit is designed for qualitative detection of antigens in fixed paraffin-embedded tissue sections, in frozen tissue sections, and in cytological samples. It is especially developed for double colour immunostaining with pairs of antibodies, one from mice, and one from rabbit.

The kit can be used for examining tissues fixed in different solutions, e.g. formalin (neutrally buffered), B5, Bouin, ethanol, or HOPE.

It is intended for research use only.

Summary and Explanation

The purpose of immunohistochemical staining is to make tissue and cell antigens visible.

ZytoChem Plus Double Stain Polymer Kit is a highly sensitive detection kit intended for use in immunohistochemistry and immunocytochemistry. The enzyme polymers in this kit consist of several molecules of secondary antibodies covalently bound to several enzyme molecules.

The Double Stain Polymer Kit includes two types of enzyme polymers: (1.) horse radish peroxidase (HRP) bound to anti-mouse secondary antibodies and (2.) alkaline phosphatase (AP) bound to anti-rabbit secondary antibodies.

Visualisation occurs via two consecutive enzyme-substrate reactions in the presence of colourising reagents which permit microscopical analysis.

The first enzyme-substrate reaction is HRP with DAB/H₂O₂. All bound primary mouse antibodies are stained in dark brown colour. The second reaction is AP with a suitable chromogenic substrate detecting all bound primary antibodies from rabbit. The chromogenic substrate for this step should result in a good contrast to DAB. We recommend using Permanent AP Red (ZUC001).

The test system is suitable for the detection of mono- and polyclonal primary antibodies and sera obtained from mice and rabbits. In contrast to other detection techniques, which often use the streptavidin-biotin technique all ZytoChem Plus Polymer Kits avoid the problem of background staining caused by endogenous biotin in the tissue.

Principle of the method

Paraffin-embedded tissue sections are first deparaffinised and rehydrated.

Endogenous peroxidase activity in the tissue may cause non-specific staining. This enzyme activity can be blocked by incubation with 3% H₂O₂-solution (peroxide block).

Background staining caused by unspecific binding of the primary antibody or the secondary antibody in the HRP-polymer is minimized by incubation with a protein blocking solution. This step can be omitted if the primary antibodies are diluted in an appropriate buffer.

The next step is the incubation with the specific primary antibodies; one must be from mice and the other from rabbit.

The antibodies can be incubated simultaneously (as a mixture or "cocktail") or one after another.

After washing, the prepared polymer mixture/working solution is applied and incubated. Any excess of unbound polymer is thoroughly washed away thereafter. The addition of DAB/H₂O₂ solution starts the enzymatic reaction of the peroxidase which leads to colour precipitation where the primary mouse antibody is bound.

The colour can be observed with a light microscope.

After another wash step the AP reaction is started by adding the second chromogenic substrate. The formed precipitate detects all bound primary rabbit antibodies.

The chromogen DAB forms a dark brown precipitate at the place of the target antigen. Permanent AP Red (compatible with AP) leads to the formation of a magenta-red product of reaction at the place of the target antigen.

Reagents provided

REF / Cat. No. POLDS-006

6 ml AP-Polymer x Rabbit (ready-to-use)
120 µl HRP-Polymer x Mouse (50 X concentrate)
1 Dilution Vial

REF / Cat. No. POLDS-100

100 ml AP-Polymer x Rabbit (ready-to-use)
2 ml HRP-Polymer x Mouse (50 X concentrate)
1 Dilution Vial

Substrate systems recommended:

DAB Substrate Kit	Cat. No. DAB057	500 tests
(2 components)	Cat. No. DAB530	5000 tests
Permanent AP Red	Cat. No. ZUC001-125	1250 tests
	Cat. No. ZUC001-500	5000 tests

Other available substrate systems:

DAB High Contrast Kit	Cat. No. DAB500 plus	500 tests
(2 Components)	Cat. No. DAB5000 plus	5000 tests
AEC Single Solution	Cat. No. ZUC037-008	80 tests
	Cat. No. ZUC037-125	1250 tests
AEC Substrate Kit	Cat. No. ZUC042-050	500 tests
(2 Components)	Cat. No. ZUC042-500	5000 tests
Permanent AEC Kit	Cat. No. ZUC054-200	2000 tests
Permanent HRP Green	Cat. No. ZUC070-100	1000 tests

Materials required but not supplied

Positive and negative control tissue
Xylene or suitable substitutes
Ethanol, distilled or deionised H₂O
Reagents for enzyme digestion or heat pre-treatment
Wash buffer PBS or TBS (Cat. No. ZUC020)
PAP Pen (Cat. No. LP0001)
Primary antibodies or antibody cocktail (user-defined)
Primary antibody diluent (Cat. No. ZUC025)
Negative control reagent
Chromogenic substrate (1 x HRP-compatible, 1 x AP-compatible)
Counterstain solution
Mounting medium
Cover slips

Storage and handling

The solutions should be stored at 2-8°C without further dilution. Please store the reagents in a dark place and do not freeze them. Under these conditions the solutions are stable up to the expiry date indicated on the label. They should not be used after the expiry date.

A positive and a negative control have to be carried out in parallel to the test material. If you observe unusual staining or other deviations from the expected results which could possibly be caused by the kit reagents, please contact Zytomed Systems' technical support or your local distributor.

Precautions

Use by qualified personnel only.

Wear protective clothing to avoid eye, skin or mucous membrane contact with the reagents. In case of a reagent coming into contact with a sensitive area, wash the area with large amounts of water. Sodium Azid and ProClin 300 are used for stabilisation. Material safety data sheets (MSDS) are available upon request. Microbial contamination of the reagents must be avoided, since otherwise non-specific staining might appear.

Reagent preparation

- 1) Reagents should be at room temperature when used.
- 2) Pipette 1 ml **AP-Polymer x Rabbit** into the provided dilution vial.
- 3) Add 20 µl **HRP-Polymer x Mouse**. Mix thoroughly.

(Note: Centrifuge the vial with the concentrated HRP-Polymer before opening.)

The prepared working solution is stable for at least 2 weeks if stored at 2–8°C. If required, smaller or larger amounts of working solution can be prepared.

Typical staining procedure

- | | |
|---|-------------------|
| 1. Blocking Solution (protein block) <i>(This step is optional.)</i> | 5 min. |
| 2. Rinse with wash buffer | 1 x 2 min. |
| 3. First primary antibody (optimally diluted) or negative control reagent | 30-60 min. |
| 4. Rinse with wash buffer | 3 x 2 min. |
| 5. Second primary antibody (optimally diluted)
<i>(Both antibodies can be incubated in one step)</i> | 30-60 min. |
| 6. Rinse with wash buffer | 3 x 2 min. |
| 7. Prepared polymer working solution | 30 min. |
| 8. Rinse with wash buffer | 3 x 2 min. |
| 9. DAB <i>(Controlling the colour intensity via light microscope is recommended.)</i> | 10 min. |
| 10. Rinse with wash buffer | 3 x 2 min. |
| 11. Permanent AP Red
<i>(Controlling the colour intensity via light microscope is recommended.)</i> | 20 min. |
| 12. Stop the reaction with distilled H ₂ O when the desired colour intensity is attained | |
| 13. Counterstaining and blueing | |
| 14. Mounting: aqueous or permanent depending on the chromogenic substrates used. | |

Quality control

We recommend carrying out a positive and a negative control with every staining run. The positive control permits the validation of appropriate processing of the sample. If the negative control has a positive result, this points to unspecific staining.

Expected results

During the reaction of the substrates with the enzyme polymers (HRP or AP) in presence of chromogens, coloured precipitates are formed at the location of the bound primary antibodies. These reactions only take place if the target antigens are existent in the tissue. The chromogens used determine the colours of the precipitates. The analysis is carried out using a light microscope.

Limitations of the procedure

Immunohistochemistry is a complex method in which histological as well as immunological detection methods are combined. Tissue processing and handling prior to immunostaining, for example variations in fixation and embedding or the inherent nature of the tissue can cause inconsistent results (Nadji and Morales, 1983).

The reagent system is especially developed for double colour staining with pairs of primary antibodies. One antibody has to be from mouse, one from rabbit. Primary antibodies from other species have to be detected via different detection systems. The same applies for two primary mouse antibodies or two primary rabbit antibodies.

It is essential to check if both antibodies are compatible with the same epitope retrieval technique.

Endogenous peroxidase or pseudoperoxidase activity may cause non-specific staining. The enzyme activity is blocked by incubation with hydrogen peroxide solution. Tissues containing Hepatitis B Surface Antigen (HBsAg) may give false positive results with HRP (horse radish peroxidase) detection systems (Omata et al, 1980).

Endogenous alkaline phosphatase activity may cause non-specific staining too. The enzyme activity can be blocked by incubation with levamisole. However, neither intestinal nor placental alkaline phosphatase can be blocked with levamisole.

Inadequate counterstaining and mounting can influence the interpretation of the results. The colour intensity of the reaction product can decrease with time, especially when exposed to light.

Zytomed Systems guarantees that the product will meet all requirements described from its shipping date until its expiry date, as long as the product is correctly stored and utilized. No additional guarantees can be given. Under no circumstances shall Zytomed System be liable for any damages arising out of the use of the reagent provided.

Troubleshooting

If you observe unusual staining or other deviations from the expected results which could possibly be caused by the reagents, please read these instructions carefully, contact Zytomed Systems' technical support or your local distributor.

No staining on an actually positive control slide:

1. Reagents were not used in the proper order.
2. Chromogenic substrate solution was too old.
3. Bleaching because chromogen and mounting medium are incompatible.
4. The antigen/epitope in the tissue was insufficiently accessible to the primary antibody. Try a pre-treatment such as heat pre-treatment or enzyme digestion. If you used a pre-treatment it should be extended.
5. Primary antibody not from mouse or rabbit, but from a different species.
6. The antigen/epitope was not stable in the fixation and/or pre-treatment procedure used. Try another fixation or pre-treatment.

Weak staining:

1. Inadequate fixation or overfixation.
2. Incomplete deparaffinisation.
3. The antigen/epitope in the tissue was insufficiently accessible to the primary antibody. If you used heat pre-treatment or enzyme digestion it should be extended.
4. Excessive incubation with Blocking Solution or insufficient washing after this step.
5. Too much wash buffer remains on the slides after washing, diluting the reagents applied in the next step.
6. If you are using PBS-based wash buffer: the activity of alkaline phosphatase in the reagents is blocked if too much wash buffer remains on the slides.
7. Incubation times were too short or primary antibody concentration too low.
8. Chromogenic substrate solution was too old.

Non-specific background staining or overstaining:

1. Incomplete deparaffinisation.
2. Excessive tissue adhesive on slides.
3. Insufficient washing especially after the incubation with the enzyme polymer mixture or the chromogenic substrate solution. These washings are critical.
4. Tissue was allowed to (partially) dry out with reagents on.
5. Unspecific binding of the primary antibodies. Please use the Blocking Solution provided with this kit or dilute the primary antibodies in appropriate diluents.
6. Incubation times of the primary antibodies were too long or primary antibody concentration too high.
7. Incubation times of the chromogenic substrate solutions were too long or reaction temperature too high (e.g. if temperature in the laboratory is high).
8. The substrate for the horse radish peroxidase is metabolised by endogenous HRP in the tissue. Maybe the hydrogen peroxide solution used for blocking was inactivated.
9. The substrate for the alkaline phosphatase is metabolised by endogenous AP in the tissue. This undesired activity can often be suppressed using levamisole (see section Limitations of the Procedure).

Performance characteristics

Zytomed Systems has conducted studies to evaluate the performance of the kit reagents. The product has been found to be suitable for the intended use.

Bibliography

Elias JM Immunohistopathology – A practical Approach to Diagnosis ASCP Press 2003
 Nadjji M and Morales AR Ann N.Y. Acad Sci 420:134-139, 1983
 Omata M et al. Am J Clin Pathol 73: 626-632, 1980

February 13, 2019

Rev: A0219

Doc: DBE_POLDS-006_100

Explanation of the symbols on the product label:

	Bestellnummer Catalog Number Reference du catalogue		Verwendbar bis Use By Utiliser jusque			Gebrauchsanweisung beachten Consult Instructions for use Consulter les instructions d'utilisation
	Chargenbezeichnung Batch Code Code du lot		Lagerungstemperatur Temperature Limitation Limites de température			Nur für Forschungszwecke For Research Use Only Pour la recherche uniquement
	In vitro Diagnostikum In Vitro Diagnostic Medical Device Dispositif médical de diagnostic in vitro		Achtung/Gefahr Warning/Danger Attention/Danger		 Hersteller / Manufacturer / Fabricant Zytomed Systems GmbH Anhaltinerstraße 16 14163 Berlin, Germany Tel: (+49) 30-804 984 990 www.zytomed-systems.de	
	Achtung/Gefahr Warning/Danger Attention/Danger		Achtung Warning Attention			